


19.09. HIS preferences

Organisation name: Enter the hospital name here.

Country: Enter the country you are located in.

Max picture size: Enter the maximum size of patient pictures that can be inserted. We suggest you set this value to 500KB.

Maximum number of users: Enter the maximum number of users that can log into the HIS system at one time. For example, if you purchase a 10 user licence, and assign 6 users to the HIS system, it will allow 4 users to log in to the pharmacy (mSupply) module.



The screenshot shows a window titled "HIS Preferences" with a standard macOS-style title bar (red, yellow, green buttons). Inside the window, there are four input fields arranged vertically on the left side:

- Organization Name: A text input field.
- Country: A text input field.
- Max picture size: A text input field containing the number "0", followed by the text "KB".
- Max number of users: A text input field.

At the bottom right of the window, there are two buttons: "Cancel" and "OK".

Previous: [19.08. Statistics](#) | | Next: [20. Integration of Other Systems with mSupply](#)

From:

<https://docs.msupply.org.nz/docs/> - **mSupply documentation wiki**

Permanent link:

<https://docs.msupply.org.nz/docs/his:preferences?rev=1584353998>

Last update: **2020/03/16 10:19**

